

LAMBDA METHAMETER

measures the concentration of CH₄ in the outgas



METHAMETER measures the **concentration of the CH₄** in outgas/exit gas.

Connected to the outgas condenser of the MINIFOR laboratory fermentor – bioreactor, the METHAMETER measures the concentration of the CH₄ gas produced during fermentation.

For more information and questions, please contact us at support@lambda-instruments.com

TECHNICAL SPECIFICATIONS

- CH₄ concentration: 0 - 9.99% or 10.0 - 99.9%
- Universal power supply: 12V/6W
- Dimension: 105 x 80 x 170 mm (W x H x D)
- Weight: 600 g
- RS communication
- CE conform